

EasyFlow™

Filtered Suction for Orthopaedic Surgery

NEW!



Unique & innovative 'plug & go' Filter

Helps to remove large bone fragments during orthopaedic surgery.

Large chamber capacity

Ideal for hip procedures and other surgeries

Prevents debris blockages

Ideal for peri-operative blood collection

Quick and Easy to remove

Cleans without disconnecting the suction line.

Spare filter

Quick for exchange

Ergonomic Handle

Ensures stability and control

Angled offset outlet

Reduces risk of tube kinking

Reversible suction probe tip

Provides an extra length for intra-femoral suction

Flexi-rib Tubing

Excellent drapability without kinking.

Gateway Medical, Inc
559 Davidson Gateway Drive
Suite 201
Davidson, NC 28036
Phone 704-655-9040
Fax 704-987-2035

Email: sales@gatewaymed.com

SMD1035



EasyFlow™

Filtered Suction for Orthopaedic Surgery

NEW!

Orthopaedic Suction Sets



- **Rapid Change Filter**
- **Ergonomic Design**
- **High Capacity**
- **Low Cost**

Product	Pennine Code	Size	Box Qty
Orthopaedic suction set - 3m Flexirib tubing - 2 curved probes - 1 straight probe - 2 filters	EFS - 1111	-7mm ID/3m length	10

Gateway Medical, Inc
559 Davidson Gateway Drive
Suite 201
Davidson, NC 28036
Phone 704-655-9040
Fax 704-987-2035

Email: sales@gatewaymed.com

